

Product Test Sheet

Product Information

Product Name	Single-wavelength Laser
Product Code	STFLH-1064-40-SM-B (LL-FCL-CW-1064-40-SMF-B)
Serial Number	26031203
Test Date	2026/04/01
Test Result	Pass

Optical Parameters

Parameter	Unit	Requirement	Test Value	Test Result
Operating Wavelength	nm	1064±2	1064.18	☒ Pass ☐ Fail
Output Power	W	1~10,adjustable	1~10,adjustable	☒ Pass ☐ Fail
Monitor Power	W	-	=1% Output Power	☒ Pass ☐ Fail
Laser Linewidth@3dB	nm	≤0.1	0.04	☒ Pass ☐ Fail
Power Peak stability (P-P)	-	No requirement	±0.22%@10W, 90min	☒ Pass ☐ Fail
Power stability (RMS)	-	No requirement	0.16%@10W, 90min	☒ Pass ☐ Fail
Output Port Isolation	dB	No requirement	35	☒ Pass ☐ Fail

Product Specifications

Operation Temperature Range	°C	-5~+50	
Operation Humidity Range	%	0~70	
Storage Temperature Range	°C	-40~+85	
Power Supply	-	AC 100~240V	USA Standard
Electrical Power Consumption	W	85	25°C
Fiber Connector	-	FC/APC	
Input/Output Fiber Type	-	Hi1060 Single Mode Fiber	Φ3mm×1m
Optical Fiber Sleeve	-	Loose tube(Φ=3mm), blue color	
Dimensions	mm	376(W)×340(D)×112(H)	Bench-Top
Remote Control Port	-	RS232-DB9	
Software Version	-	LaserController-dBm-V20250222	

SIMTRUM China
Telephone: +86 133 2643 0008
Email: info@simtrum.com

SIMTRUM Singapore
Telephone: +65 6996 0391
Email: info@simtrum.com



8 Evaluation

8.1 General

The standard deviation, s , from n readings m_i is calculated according to

$$s = \sqrt{\frac{\sum_{i=1}^n (m_i - \bar{m})^2}{n-1}} \quad (3)$$

Where the mean value is


$$\bar{m} = \frac{\sum_{i=1}^n m_i}{n} \quad (4)$$

8.3 Power stability of cw lasers

Calculate the mean value of the power, $\bar{P}$, and the respective standard deviation, s , for the appropriate stability time domain (short-term, medium-short-term, medium-term and long-term) according to the specifications given in 7.3.

Power stability is given as the relative power fluctuation, ΔP , in the corresponding stability time domain calculated from [Formula \(9\)](#):


$$\Delta P = \frac{2s}{p} \quad (9)$$

Spectral Test	
Test Item	The spectrum of the Laser
Test Method	After 20dB attenuation,access to the spectrometer for measurement
Measurement Equipment	YOKOGAWA AQ6370D
Measurement Results	 YOKOGAWA // AQ6370D OPTICAL SPECTRUM ANALYZER // 2026 Mar 31 13:50 TR A VPK : 1064.1800nm 17.98dBm V-Vn: V0001: /BLK /BLK V0002: /BLK /BLK V0003: /BLK /BLK V0004: /BLK /BLK V0005: /BLK /BLK (MEAS CONDITION) START: 1059.000nm STOP: 1069.000nm CENTER: 1064.000nm SPAN: 50.0nm 10.0 dB/D RES: 0.1nm SENS: 110 AVG: 1 SML: 2501 (AUTO) REF. 0.0 dBm 1059.000nm 1064.000nm 1069.000nm 5.00 nm/D 20.0 0.0 -20.0 -40.0 -60.0 -80.0 200 100 50 20 10 5 2 1 0.5 0.2 0.1 0.05 0.02 0.01 200 100 50 20 10 5 2 1 0.5 0.2 0.1 0.05 0.02 0.01 RPT SGL STP
	Optical Spectrum@10W, Monitor Port

SIMTRUM China
Telephone: +86 133 2643 0008
Email: info@simtrum.com

SIMTRUM Singapore
Telephone: +65 6996 0391
Email: info@simtrum.com





Power Stability Test	
Test Item	Laser Power and Power Stability
Test Method	Direct Measurement with Optical Power Meter
Measurement Equipment	Thorlabs PM100D/S146C Power Meter
Measurement Results	P (W) T (min) P-P: $\pm 0.22\%$ RMS: 0.16%

Power Stability Test@10W

SIMTRUM China
Telephone: +86 133 2643 0008
Email: info@simtrum.com

SIMTRUM Singapore
Telephone: +65 6996 0391
Email: info@simtrum.com



Inspection of the fiber optic connector end face	
Test item	Integrity and cleanliness of optical fiber connector face
Test method	The optical fiber end face mirror is directly observed with a magnification of 400x.
Equipment	Hand-held End Face Detector/DL-AutoGet-V2
Test result	The output fiber end face is intact, free from scratches, and clean without dust or contamination, meeting all requirements. Output Port Monitor Port

SIMTRUM China
Telephone: +86 133 2643 0008
Email: info@simtrum.com

SIMTRUM Singapore
Telephone: +65 6996 0391
Email: info@simtrum.com

